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Caution: The use of the device is limited to licensed physicians or practitioners in medical institutions formally approved by the government with a professional training certificate issued by the manufacturer or its entrusted/ designated institutions in strict accordance with the Instructions for Use.

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

COMPOSITION

Hyaluronic Acid gel	24 mg
Phosphate buffer pH 7.2 q.s.	1 mL
One syringe contains 0.8 mL of Juvéderm ULTRA PLUS®	
*including 5 % of uncrosslinked Hyaluronic Acid.	

1. DEVICE DESCRIPTION

Juvéderm ULTRA PLUS® injectable gel is a sterile, biodegradable, nonpyrogenic, viscoelastic, clear, colorless, homogeneous gel implant. **Juvéderm ULTRA PLUS®** injectable gel consists of cross linked hyaluronic acid (HA) produced by *Streptococcus* bacteria. The gel is presented in a graduated pre-filled disposable syringe. Each box contains two 0.8 mL **Juvéderm ULTRA PLUS®** syringes, 4 single-use 27G1/2" sterile needles to be used only for injecting **Juvéderm ULTRA PLUS®**, an instruction leaflet and a set of labels in order to ensure traceability. The content of **Juvéderm ULTRA PLUS®** syringes is sterilized by moist heat. The 27G1/2" needles are sterilized by radiation.

2. INTENDED USE/INDICATIONS

Juvéderm ULTRA PLUS® injectable gel is indicated for injection into the mid to deep dermis for correction of severe nasolabial folds.

3. CONTRAINDICATIONS

- **Juvéderm ULTRA PLUS®** injectable gel is contraindicated for patients with severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies.
- **Juvéderm ULTRA PLUS®** injectable gel contains trace amounts of gram-positive bacterial proteins and is contraindicated for patients with a history of allergies to such material.

4. WARNINGS

- **Juvéderm ULTRA PLUS®** injectable gel must not be injected into blood vessels. Introduction of **Juvéderm ULTRA PLUS®** injectable gel into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting soft tissue fillers, for example, after insertion of the needle, and just before injection, the plunger rod can be withdrawn slightly to

aspirate and verify the needle is not intravascular, inject the product slowly and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate physician should an intravascular injection occur. (See Physician Instructions #10)

- Use of **Juvéderm ULTRA PLUS**[®] injectable gel at specific sites in which an active inflammatory process (skin eruptions such as cysts, pimples, rashes, or hives) or infection is present should be deferred until the underlying process has been controlled.
- Injection procedure reaction to **Juvéderm ULTRA PLUS**[®] injectable gel has been observed as consisting mainly of short-term inflammatory symptoms starting early after treatment and with less than 7 days' duration. Refer to the SIDE EFFECTS section for details.

5. PRECAUTIONS

- **Juvéderm ULTRA PLUS**[®] injectable gel is packaged for single patient use. Do not resterilize. Do not use if package is opened or damaged.
- In order to minimize the risks of potential complications, this product should only be used by physicians who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Physicians are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Based on preclinical studies, patients should be limited to 20 mL of any Juvéderm[®] injectable gel per 60 kg (130 lbs) body mass per year. The safety of injecting greater amounts has not been established.
- The safety and effectiveness of **Juvéderm ULTRA PLUS**[®] injectable gel for the treatment of anatomic regions other than severe nasolabial folds have not been established in controlled clinical studies.
- As with all transcutaneous procedures, **Juvéderm ULTRA PLUS**[®] injectable gel implantation carries a risk of infection. Standard precautions associated with injectable materials should be followed.
- **Juvéderm ULTRA PLUS**[®] is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity, and performance of the product, and it can therefore no longer be assured.
- The safety of **Juvéderm ULTRA PLUS**[®] injectable gel for use during pregnancy, in breastfeeding females, or in

patients under 18 years has not been established.

- The safety of **Juvéderm ULTRA PLUS**[®] injectable gel in patients with known susceptibility to keloid formation, hypertrophic scarring, and pigmentation disorders has not been studied.

- **Juvéderm ULTRA PLUS**[®] injectable gel should be used with caution in patients on immunosuppressive therapy.

- Patients who are using substances that can prolong bleeding (such as aspirin, nonsteroidal anti-inflammatory drugs, and warfarin) may, as with any injection, experience increased bruising or bleeding at injection sites.

- After use, treatment syringes and needles may be potential biohazards. Handle and dispose of these items in accordance with accepted medical practice and applicable local, state, and federal requirements.

- **Juvéderm ULTRA PLUS**[®] injectable gel is a clear, colorless gel without visible particulates. In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe; notify Allergan Product Support immediately.

- If laser treatment, chemical peeling, or any other procedure based on active dermal response is considered after treatment with **Juvéderm ULTRA PLUS**[®] injectable gel, there is a possible risk of eliciting an inflammatory reaction at the implant site. This also applies if **Juvéderm ULTRA PLUS**[®] injectable gel is administered before the skin has healed completely after such a procedure.

- Failure to comply with the needle attachment instructions could result in needle disengagement and/or product leakage at the LUER-LOK[®] and needle hub connection.

- In the China registration trial, the primary efficacy endpoint was measured for up to 6 months; the majority of Chinese subjects showed continued clinically meaningful improvement over baseline at 6 months.

- It is recommended that the maximum injection volume of **Juvéderm ULTRA PLUS**[®] injectable gel is 1.5 mL per NLF for initial injection and 0.5 mL per NLF for touch-up injection (4 weeks later after initial injection) which is needed for optimal correction based on China clinical study.

6. SIDE EFFECTS

- Injection procedure reactions to **Juvéderm ULTRA PLUS**[®] injectable gel has been observed as consisting mainly of short-term inflammatory symptoms such as swelling, induration, pain, bumps, redness, bruising and itching.

Subjects mainly had mild to moderate injection-site responses. All injection site reactions completely resolved without medical intervention.

7. INSTRUCTIONS FOR USE

A. To Attach Needle to Syringe

STEP 1: Remove tip cap

Hold syringe and pull tip cap off the syringe as shown in Figure A.

Figure A


STEP 2: Insert needle

Hold the syringe body and firmly insert the hub of the needle (provided in the **Juvéderm**[®] package) into the LUER-LOK[®] end of the syringe.

STEP 3: Tighten the needle


Tighten the needle by turning it firmly in a clockwise direction (see Figure B) until it is seated in the proper position, as shown in Figure C.

NOTE: If the position of the needle cap is as shown in Figure D, it is not attached correctly. Continue to tighten until the needle is seated in the proper position.

STEP 4: Remove the needle cap

Hold the syringe body in one hand and the needle cap in the other. Without twisting, pull in opposite directions to remove the needle cap as shown in Figure E.

B. Physician Instructions

1. **Juvéderm ULTRA PLUS**[®] injectable gel is a highly crosslinked formulation that can be injected using a 27-G needle for more versatility in contouring and volumizing of severe nasolabial folds. Prior to treatment with **Juvéderm ULTRA PLUS**[®] injectable gel, the patient's medical history should be obtained and the patient should be fully apprised of the indications, contraindications, warnings, precautions, treatment responses, adverse reactions, and method of administration. Patients also should be advised that supplemental "touch-up" implantations may be required to achieve and maintain maximum correction.
2. The patient's soft-tissue deficiencies should be characterized with regard to etiology, distensibility, stress at the site, and depth of lesion. Depending on the type of skin, best results are obtained when the defect is readily distensible and correction can be visualized by manual manipulation (stretching) of the skin. Pretreatment photographs are recommended.

3. Topical or injectable anesthesia may be used to manage pain during and after injection.
4. After ensuring that the patient has thoroughly washed the treatment area with soap and water, the area should be swabbed with alcohol or other antiseptic. Prior to injecting **Juvéderm ULTRA PLUS**[®] injectable gel, depress the plunger rod until the product flows out of the needle.
5. The injection technique of **Juvéderm ULTRA PLUS**[®] injectable gel with regard to the angle and orientation of the bevel, the depth of injection, and the quantity administered may vary. A linear threading technique, serial puncture injections, or a combination of the two have been used to achieve optimal results. If **Juvéderm ULTRA PLUS**[®] injectable gel is injected too superficially, this may result in visible lumps and/or discoloration.
6. Inject **Juvéderm ULTRA PLUS**[®] injectable gel applying even pressure on the plunger rod while slowly pulling the needle backward. The wrinkle should be lifted and eliminated by the end of the injection. It is important that the injection is stopped just before the needle is pulled out of the skin to prevent material from leaking out or ending up too superficially in the skin.
7. If the needle is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle.
8. The typical total volume of **Juvéderm ULTRA PLUS**[®] injectable gel to achieve optimal correction for initial treatment is limited to 1.5 mL per NLF. The typical volume of **Juvéderm ULTRA PLUS**[®] injectable gel to achieve optimal correction for touch-up injection is limited to 0.5 mL per NLF.
9. Correct to 100% of the desired volume effect. Do not overcorrect. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue, and the injection technique. Markedly indurated defects may be difficult to correct
10. If immediate blanching occurs, the injection should be stopped and the area massaged until it returns to a normal color. Blanching may represent a vessel occlusion. If normal skin coloring does not return, do not continue with the injection.
11. When injection is completed, the treated site should be gently massaged so that it conforms to the contour of the surrounding tissues. If overcorrection occurs, massage the area between your fingers or against an underlying superficial bone to obtain optimal results.
12. With patients who have localized swelling, the degree of correction is sometimes difficult to judge at the time of treatment. In these cases, it is better to invite the patient to a touch-up session after 4 weeks.
13. Patients may have mild to moderate injection-site responses, which typically resolve in a few days. If the treated area is swollen immediately after the injection, an ice pack can be applied to the site for a short period.
14. After the initial treatment, an additional treatment of **Juvéderm ULTRA PLUS**[®] injectable gel (from 4 weeks later) may be necessary to achieve the desired level of

correction. If the wrinkle needs further treatment, the same procedure should be repeated until a satisfactory result is obtained. The need for an additional treatment may vary from patient to patient and is dependent upon a variety of factors such as wrinkle severity, skin elasticity, and dermal thickness at the treatment site.

15. The physician should instruct the patient to promptly report to her/him any evidence of problems possibly associated with the use of **Juvéderm ULTRA PLUS®** injectable gel.

C. Patient Instructions

It is recommended that the following information be shared with patients:

- Within the first 24 hours, patients should avoid strenuous exercise, extensive sun or heat exposure, and alcoholic beverages. Exposure to any of the above may cause temporary redness, swelling, and/or itching at the injection sites.
- To report an adverse reaction, contact the Allergan Product Support.

8. HOW SUPPLIED

Juvéderm ULTRA PLUS® injectable gel is supplied in individual treatment syringes with 27-G needles for single-patient use and ready for injection (implantation). The volume in each syringe is as stated on the syringe label and on the carton. The contents of the syringe are sterile and nonpyrogenic. Do not resterilize. Do not use if package is opened or damaged.

9. SHELF LIFE AND STORAGE

Juvéderm ULTRA PLUS® injectable gel must be used prior to the expiration date printed on the label.

Store at room temperature (up to 25°C/77°F). DO NOT FREEZE.

Juvéderm ULTRA PLUS® injectable gel has a clear appearance. In the event that a syringe contains material that is not clear, do not use the syringe; notify Allergan Product Support immediately.

To place an order, contact Allergan.

Product shelf life: 2 years

Date of manufacture : refer to original label.

Product name:	Juvéderm ULTRA PLUS®
Product Model and specifications :	0.8ml/Syringe 2 Syringes / box
Registration certificate number:	GXRI 20153131708
Product technical requirements :	GXRI 20153131708
Manufacturer / Registrant name:	Allergan
Manufacturer / Registrant address:	Route de Promery Zone Artisanale de Pré-Mairy Pringy 74370 Annecy FRANCE
Manufacturing Site address:	Route de Promery Zone Artisanale de Pré-Mairy Pringy 74370 Annecy FRANCE
Manufacturer / Registrant Tel:	+ 33(0) 4 50 27 27 03
Allergan Product Support e-mail:	MedDeviceComplaintsAPAC@Allergan.com
Agent name:	ALLERGAN INFORMATION CONSULTING (SHANGHAI) CO.,LTD
Agent Address:	Room 5605, 1266 West Nanjing Road, Jingan District, Shanghai, China 200040
Agent Tel:	021-60329100
Agent Fax:	021-60329104
After-sales service facility:	SINOPHARM GROUP DISTRIBUTION CO.,LTD.
After-sales service Address:	No.270 Mei Yue Road China (Shanghai) Pilot Free Trade Zone 200131
After-sales service Tel:	4001538070
After-sales service Fax:	021-51062770

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